

## Product Datasheet

### Goat anti-dengue antibody #SB-40-3020

<b>Immunogen</b>	10- amino acid long internal peptide from the Dengue virus
<b>Quantity</b>	0.2mg in 0.5mg/ml
<b>Purification</b>	Ammonium sulphate precipitation of serum followed by antigen affinity purification
<b>Buffer</b>	Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative
<b>Storage and handling</b>	Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles
<b>Applications tested</b>	Peptide ELISA: Antibody detection limit dilution 1:64000
<b>Application notes</b>	Please optimize dilutions for your particular assay conditions.