

Product Datasheet

Goat anti-dengue antibody #SB-40-3021

Immunogen	10- amino acid long internal peptide from the Dengue virus
Quantity	0.2mg in 0.5mg/ml
Purification	Ammonium sulphate precipitation of serum followed by antigen affinity purification
Buffer	Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative
Storage and handling	Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles
Applications tested	Peptide ELISA: Antibody detection limit dilution 1:64000
Application notes	Please optimize dilutions for your particular assay conditions.