

## Product Datasheet

### **Goat anti-Human HCG antibody #SB-20-3001**

<b>Immunogen</b>	a 35-aa long internal peptide from the HCG hormone
<b>Quantity</b>	0.2mg in 0.5mg/ml
<b>Purification</b>	Ammonium sulphate precipitation of serum followed by antigen affinity purification
<b>Buffer</b>	Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative
<b>Storage and handling</b>	Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles
<b>Applications tested</b>	Peptide ELISA: Antibody detection limit dilution
<b>Application notes</b>	Please optimize dilutions for your particular assay conditions.