

Product Datasheet

Goat anti-Human HCG antibody #SB-20-3002

Immunogen	a 35-aa long internal peptide from the HCG hormone
Quantity	0.2mg in 0.5mg/ml
Purification	Ammonium sulphate precipitation of serum followed by antigen affinity purification
Buffer	Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative
Storage and handling	Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles
Applications tested	Peptide ELISA: Antibody detection limit dilution
Application notes	Please optimize dilutions for your particular assay conditions.