

Product Datasheet

Goat anti-dengue antibody #SB-40-3020

| | |
|-----------------------------|--|
| Immunogen | 10- amino acid long internal peptide from the Dengue virus |
| Quantity | 100µg in 200µl |
| Purification | Ammonium sulphate precipitation of serum followed by antigen affinity purification |
| Buffer | Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative |
| Storage and handling | Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles |
| Applications tested | Peptide ELISA: Antibody detection limit dilution 1:64000 |
| Application notes | Please optimize dilutions for your particular assay conditions. |