

Product Datasheet
Goat anti-Human HCG antibody
#SB-20-3001

Immunogen a 35-aa long internal peptide from the HCG hormone

Quantity 100µg in 200µl

Purification Ammonium sulphate precipitation of serum followed by antigen affinity purification

Buffer Phosphate buffer saline, PH 7.2, with 0.02% sodium azide as preservative

Storage and handling Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles

Applications tested Peptide ELISA: Antibody detection limit dilution

Application notes Please optimize dilutions for your particular assay conditions.