

Product Datasheet

Goat anti-dengue antibody #SB-40-3021

Immunogen	10- amino acid long internal peptide from the Dengue virus
Quantity	100µg in 200µl
Purification	Ammonium sulphate precipitation of serum followed by antigen affinity purification
Buffer	Phosphate buffer saline, PH 7.2, with 0.5% BSA (US origin) and 0.02% sodium azide
Storage and handling	Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles
Applications tested	Peptide ELISA: Antibody detection limit dilution 1:64000
Application notes	Please optimize dilutions for your particular assay conditions.